

# SLIDEVIEW™ DX VS-M1

## Whole Slide Imaging System

### Fast, Efficient Digital Pathology



**EVIDENT**

Our digital pathology slide scanner is sold as a research use only device in the United States and China and as a medical device in Europe, Singapore, Australia, New Zealand, India, Korea, and Japan. Product currently not available in regions not specified.

# High-Quality Images at High Speed

We combined our leading-edge technologies—including X Line™ objectives and True Color LED illumination—to create a fast slide scanner that delivers microscope-quality images onscreen. The SLIDEVIEW™ DX whole slide imaging system's advanced technology improves the entire digital pathology workflow to help pathologists make diagnoses efficiently with high-quality images they can rely on.



This product is sold as a research use only device in the U.S. and China and should not be used for diagnostic purposes in these regions.





# Supports Efficient Diagnosis

Make your diagnoses with the convenience of digital slide images that you can have confidence in. The SLIDEVIEW™ DX system uses a dedicated objective lens based on our award-winning X Line™ series to capture images with exceptional flatness and high resolution, enabling seamless stitching and stunning whole slide images that are fully in focus.



## Microscope-Quality Images Onscreen

The objective is paired with our high color rendering True Color LED light source to provide the color reproduction necessary for H&E, IHC, and special stains.

The images are displayed on a high color rendering monitor that has sRGB wide gamut coverage tailored for digital pathology images and that supports the correctness and stability of colors over time. You can be confident that the slide images on the screen are as close as possible to what you would see looking through the oculars.



This product is sold as a research use only device in the U.S. and China and should not be used for diagnostic purposes in these regions.

## Fast, Efficient Pathology Diagnosis Workflow

The SLIDEVIEW™ DX software was developed in collaboration with pathologists, so its layout and functionality will be intuitive and easy to learn.

The system's software integrates with your lab information system (LIS) to display the slide and patient information in one convenient window to help you remain focused on making a diagnosis. A virtual slide tray mirrors handling a physical slide on the holder for simple, familiar slide management. View images from the virtual slide tray and easily make measurements and annotations on the digital slide image.





Image management system with virtual slide tray display



Annotation function

For consultation, the system has convenient collaboration tools that enable you to securely share slide images with other pathologists or experts. Users can access the image management system from a web browser without installing software.



Convenient collaboration functions



# Built for Busy Labs

The SLIDEVIEW™ DX system makes it easy to scan large numbers of slides. Labs can set up their system based on their preferences. You can configure the system for greater speed and slide quality, giving you flexibility over the scan plan and making it easy to find the right balance of speed and quality. If you need more capacity to run slides, additional scanners can be linked together. For example, using two scanners would allow you to process 300 slides at a time.

## Easy to Operate

The scanner's controls are easily accessible from its large touch screen, which is designed to be used while standing. This makes it easy for you to configure the scanner's controls and walk away.

Easily load the slides into the rack and then place them into the scanner. The slides are loaded vertically rather than horizontally to reduce the chance of dropping them. The system walks you through the steps to start scanning. When the scan is complete, you can check the quality of the digital slides using the same touch screen and deliver them directly to the pathologist at the push of a button.




Easy slide loading



Simple touch screen operation

## High-Speed Scanning


The scanner's fast, more than 80 slides per hour scan speed enables slides to be quickly digitized for faster delivery to the pathologist. During the scan, the system's real-time autofocus enables you to skip the time-consuming focus mapping phase for faster scans.



## Reliable

The system's advanced AI\* automatically identifies the tissue on the slide so that only the tissue is scanned. This speeds up the scanning process and minimizes the need for rescans caused by undetected tissue.


\* Machine learning-based software with regulatory framework.



Properly detected samples

## Flexible

The system has three scan modes that enable you to tailor the scan based on your lab's priorities. Speed scan mode uses real-time autofocus to acquire quality images quickly. The focus scan mode uses single-point autofocus to create a detailed focus map before scanning, improving accuracy. Lastly, with manual setting mode, you have the option to individually set each focus point, helping make faintly stained or uneven samples completely in focus.



In speed scan mode, the scanner automatically focuses on each grid square and captures images one by one. The images are then stitched together to create one large, all-in-focus image.

## Easy to Integrate and Manage


The Evident Management Portal (EMP) is used to configure and maintain our scanners. The EMP can be connected to multiple scanners in the laboratory. It enables you to control who has access to the scanner and EMP, set scanner configurations, such as barcode parsing and the default scan mode, and provides diagnostic information, including active status, instrument logs, audit logs, and scanner usage statistics.

The software easily integrates with your existing lab information system (LIS) and supports various healthcare network protocols, including HL7. The DICOM\* format enables compatibility with existing PACS- or DICOM-based image management systems. The system has HIPAA and NIST data security features, encryption, and central account management.

IT technicians can use the system's central account manager to remotely view scanner system faults and errors, as well as download reports.

\*Optional software solution.

## SLIDEVIEW DX Dimensions



unit: mm

## SLIDEVIEW DX Specifications

Intended Specimen	Specimen	Glass slide with cover glass
	Size of Glass Slide	Length: 75–76.5 mm (2.954–3.011 in.), Width: 25–26.5 mm (0.985–1.043 in.), Thickness: 0.9–1.2 mm (0.0355–0.0472 in.) (ISO 8037-1)
	Thickness of Cover Glass	0.13–0.19 mm (0.0052–0.0074 in.) (ISO 8255-1) (recommendation: 0.17 mm)
Scanner	Touch Screen	Built-in 21.5-inch color LCD Adjust the scan settings and observe acquired images using the touch screen
	Slide Capacity	Up to 150 slides
	Slide Rack	5 slide racks, up to 30 slides per rack
	Slide Rack Adapter	For Sakura Finetek Japan Co., Ltd. slide basket (K1650009), transferring 10 slides to the slide rack at once
	Scan Modes	Speed scan mode, focus scan mode, and manual setting mode
	Throughput	More than 80 slides/hour (15 mm × 15 mm area at 40X equivalent resolution, speed scan mode)
	Objective Lens	X Line extended apochromat
	Resolution	0.23 μm/pixel (equivalent to 40X)
	Camera	Image sensor: 1.1 inch CMOS, pixel pitch: 3.45 × 3.45 μm
	Illuminator	High intensity and high color rendering LED (up to 50,000 hours)
	Barcode Reader	1D Barcode WPC (JAN/EAN/UPC-A/UPC-E), NW-7, ITF, Industrial2of5, Code39, Code128, RSS-14, RSSLimited, RSSExpanded
		2D Barcode QR Code, DataMatrix (ECC200), MicroQR, PDF417, MicroPDF417
	Automatic Slide Detection	Yes, the slide rack and slide positions are automatically detected and displayed on the touch screen
	Automatic Scan Area Detection	Yes
	Anti-Vibration Function	Yes
	LED Status Display	Hardware status, slide rack status, scanning status
	External Dimensions (W × D × H)	1060 × 625 × 710 mm (41.7 × 24.6 × 27.9 in.)
	Weight	150 kg (330.7 lb)
	Primary Power Source	100–240 VAC, 50-60 Hz
	Power Consumption	Max. 120 W
Image Management System	Environmental Conditions	Temperature range: A maximum operating range of 7°C (44.6 °F) within the limits of 15–27 °C (59–80.6 °F) (no sudden temperature changes) Humidity range: 35–80% RH (non-condensing indoor use) Pressure range: sea level to 2000 m (6561.7 ft)
	Display Functions	Search filters, slide list, virtual slide tray, full-screen, split, macro image, slide properties
	View Functions	Pan, zoom, rotate, synchronization, slide navigator, rotation compass, heat map, scale, magnification stop, magnifying loupe
	Annotation and Measurement	Flag, square, circle, line, comment description
Software Options	Collaboration Functions	Chat, shared view
	DICOM	Scans are produced according to DICOM image standards for integration with PACS- and DICOM-based image management systems
	LIS Connector	Integration with LIS systems to present samples and patient data in the Evident image management system

• In Europe, SLIDEVIEW DX VS-M1 is regulatory compliant (under the name VS-M1) for IVDR (EU) and MDR(UK) for in vitro diagnostic use.  
 • In Singapore, SLIDEVIEW DX VS-M1 is listed (under the name of VS-M1) as a Class A medical device in HSA Class A SMDR.  
 • In Australia, SLIDEVIEW DX VS-M1 is included (under the name VS-M1) in the ARTG as a Class 1 IVD medical device.  
 • In New Zealand, SLIDEVIEW DX VS-M1 is regulatory compliant (under the name VS-M1) for the Medicines Act of 1981 (NZ).  
 • In India, SLIDEVIEW DX VS-M1 is regulatory compliant (under the name VS-M1-IVD) as a Class A medical device in CDSCO.  
 • In Korea, SLIDEVIEW DX VS-M1 is regulatory compliant (under the name VS-M1-IVD) as a Class 1 medical device in MFDS.

• **EVIDENT CORPORATION is ISO14001 certified.**  
 For details on certification registration, visit [olympus-lifescience.com/en/support/iso/](https://olympus-lifescience.com/en/support/iso/)  
 • **EVIDENT CORPORATION is ISO9001 certified.**  
 • **CLASS 1 Laser Product**  
 • **All company and product names are registered trademarks and/or trademarks of their respective owners.**  
 • **Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.**

This product is sold under the name SLIDEVIEW VS-M1R in the U.S. and China as a research use only product and as SLIDEVIEW VS-M1J in Japan notified (under the name VS-M1-IVD1) as a Class 1 medical device. Product currently not available in regions not specified.

**EvidentScientific.com**

**EVIDENT**

**EVIDENT CORPORATION**  
 Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0910, Japan

**OLYMPUS**