

Compact
Powerful
Precise



ColorFlex[®] L2 Tomato

*Level UP your Color Confidence with the new
ColorFlex L2 Tomato Spectrophotometer.*

As simple as you want, as powerful as you need.

Elevate your color accuracy and user experience like never before. Compact and ergonomic build requires minimal space while offering enhanced durability and reliability. Designed to streamline and expedite your quality control processes, ColorFlex L2 Tomato Spectrophotometer is expertly crafted to measure the color of tomatoes in all their processed forms, including paste/purée, sauce, ketchup, juice, and fresh tomatoes.




Explore the Next-Gen ColorFlex®...L2 Tomato



#1 Color Measurement Solution for Tomatoes

HunterLab Tomato Reference Standard: HunterLab, the creators of the original tomato color reference standard, offers the only reference tile that traces back to the discontinued BCR and USDA standards, ensuring the most accurate red color measurements.




All the Tomato Scores and Indices You Need

Measure Tomato in Every Form: Tomato Paste Score (TPS), Tomato Sauce Score (TSS), Tomato Catsup Score (TCS), Tomato Juice Score (TJS), a/b ratio, Fresh Tomato Color Index (FTCI) & Lycopene index.

45°a:0° Annular Illuminated Optical Geometry

Superior Consistency and Color Accuracy: Enhancing the efficiency of Quality Control workflows and delivery of superior quality products repeatedly.



Compact with Sealed Spillproof Case

Space-Saving and Easy Maintenance: Occupies minimal space and is designed to withstand the demands of busy lab environments.



Tomato Paste
Tomatoes

Tomato Sauce
TOMATO KETCHUP

tomato juice
Tomato Soup



Powerful Standalone Functionality

No computer required: With huge onboard processing and storage capacity, ColorFlex L2 Tomato can function as an independent color workstation. Connect monitor, keyboard, and mouse if desired.




Modern Tablet User Interface with Touchscreen Color Display

Durable and Efficient: Facilitates accurate color measurement in seconds, with a robust industrial glass design suitable for harsh environments.




New & Improved Color Software

Essentials 2.0: Take precise and reliable color measurement in seconds with next-gen user-friendliness and seamlessly analyze multiple tristimulus and spectral data sets. It's that powerful.



Self-Guided and Intuitive Out-of-Box Startup

Time-Saving and User-Friendly: Ready to use straight out of the box with an integrated wizard for quick training, reducing operator time and easing concerns related to staff training and turnover.



Built-in Camera for Sample Viewing and Screen Recording

Precision and Communication: Ensures accurate sample positioning to eliminate measurement discrepancies and allows easy sharing of recordings for swift quality control decisions.



Tomato Paste
Tomatoes


Tomato Sauce
TOMATO KETCHUP

tomato juice
Tomato Soup



Built-in Advanced Diagnostics

On-Board Diagnostic Health Check: Ensures the HunterLab ColorFlex L2 Tomato is always operating at peak performance, continuing to build on the legendary history of reliability of the ColorFlex family of instruments.



Advanced Data Export and Communication

Full suite of data communication options: Features USB ports, HDMI, and support for mouse and keyboard, enabling easy data export/printing.

Connectivity with LIMS and SPC systems

Simplifies QC workflows: Aids in trend analysis and workflow simplification across manufacturing sites.



CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, request a quote or just say hello!

OUR OFFICES & DISTRIBUTION CENTERS

HunterLab
 11491 Sunset Hills Road
 Reston, VA 20190-5280
 Website: www.hunterlab.com
 Email: sales@hunterlab.com
 Ph: 703-471-6870
 Fax: 703-471-4237

Corporate Headquarters

HunterLab Europe
 Christian Jansen, Managing Director
 Dr.-August-Einsele-Ring 15
 82418 Murnau, Germany
 Website: www.hunterlab.com
 Email: Christian.Jansen@hunterlab.com
 Ph: +49 8841 9464

European Office

HunterLab Asia-Pacific
 Heng Su Kiang, APAC Sales Manager
 28-2, Plaza KLTS
 No. 99, Jalan Gombak
 Setapak
 53000 Wilayah Persekutuan Kuala Lumpur
 Malaysia
 Website: www.hunterlab.com
 Email: sales-asiapacific@hunterlab.com
 Ph: +(60) 340 324 749

HunterLab Asia-Pacific
 Leo Cheung, Regional Manager
 Flat E, 13F Leahander Centre
 28 Wang Wo Tsai Street Tsuen Wan
 New Territories Hong Kong
 Website: www.hunterlab.com
 Email: sales-asiapacific@hunterlab.com
 Ph: +852-2406 1723

Asia-Pacific Offices

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com

SPECIFICATIONS



ColorFlex[®] L2 Tomato

MEASUREMENT

Measurement Principle:	Port up dual-beam spectrophotometer
Geometry:	45°a:0° Annular Illuminated Optics, ASTM E1164
Spectrophotometer:	Sealed optics; 256-element diode array and a high resolution concave holographic grating
Port Diameter/View Diameter:	31.8 mm (1.25 in) illuminated/ 25.4 mm (1 in) measured
Specular Component:	Excluded
Spectral Range:	Range: 400 nm - 700 nm
Spectral Resolution:	< 3 nm
Effective Bandwidth:	10 nm equivalent triangular
Reporting Interval:	10 nm
Photometric Range:	0 to 150 %
Light Source:	Full spectrum Xenon Flash Lamp
Measurement Duration:	1 second
Lamp Life:	10 years typical
Illumination:	Factory Calibrated Xenon for Brightness, Whiteness, and Tint
Standards Conformance:	CIE 15:2018, ASTM E1164, DIN 5033, Teil 7 and JIS Z 8722 Condition C
Standards Traceability:	Instrument standard assignment in accordance with National Institute of Standards and Technology (NIST) following practices described in CIE Publication 44 and ASTM E259

PERFORMANCE

Inter-Instrument Agreement:	Color: $\Delta E^*_{2000} < 0.15$ D65/10 (Avg) on BCRA CCSII Tile
Colorimetric Repeatability:	Color: $\Delta E^* < 0.05$ D65/10 (Max) on White Tile



USER INTERFACE

Data Views:	Color Data, Spectral Plot, EZ View, Tristimulus Color Plot, Pass/Fail Color Indication, Time, and Date Stamp, Auto-naming, Auto-saving, Data backup, Image Capture via High-resolution Camera
USB Flash Drive Features:	Job Data, Sample Images, Screen Capture and Database Export
Other Features:	Pass/Fail, Average Multiple Measurements, Set Multiple Automatic Measurements
Illuminants:	A, C, D50, D55, D65, D75, F02, F07, F11
Observers:	2° and 10°
Color Scales:	CIE L*a*b*, Hunter Lab, CIE L*C*h, CIE Yxy, CIE XYZ, Rd,a,b
Color Differences:	$\Delta L^*a^*b^*$, ΔLab , ΔL^*C^*h , ΔYxy , ΔXYZ , ΔE , ΔE^* , ΔE^*2000 , ΔH^* , ΔC^* , ΔE_{cmc}
Indices and Metrics:	E313 Whiteness and Tint (C/2°, C/10°, D65/2° or D65/10°), E313 Yellowness (C/2°, C/10°, D65/2° or D65/10°), D1925 Yellowness (C/2°), Y Brightness, Z%, 457 nm Brightness, Opacity, Color Strength (Average and Single Wavelength), Gray Scale, Gray Stain, Metamerism Index
Tomato Scores and Indices:	Tomato Paste Score (TPS) Tomato Juice Score (TJS) Tomato Sauce Score (TSS) Tomato Catsup Score (TCS) Tomato Dice Score (TDS) Tomato a/b Ratio (a/b) Lycopene Index (Lyc) Fresh Tomato Color Index (FTCI)
Data Storage:	1 million Records max; 8 GB
Languages:	Supports multiple languages—with the option to add more as needed

PHYSICAL / ELECTRICAL

Dimensions:	Height: 16 cm (6.25 in.) Width: 21 cm (8.75 in.) Depth: 24 cm (9.50 in.) Weight: 2.75 kg (6 lbs)
Display:	7" Touch screen, high resolution 1280x800
Interface:	2 USB ports, HDMI Video Output, Foot switch Input, Ethernet connection, USB Service Port, Convenient physical measurement/action button
Power:	90 – 240 VAC, 47 – 63 Hz to universal power supply @ 24 VDC/3.75A
Operating Environment:	10° to 30°C (50° to 86° F), 10 % to 90 % RH, noncondensing
Storage Environment:	-20° to 65°C (-4° to 149° F), 10 % to 90 % RH, noncondensing
Standard Accessories:	<ul style="list-style-type: none"> • HunterLab Reference Standard-Certified Tomato with Certificate of Traceability • Calibrated Instrument Standard (White Tile) with Certificate of Traceability • Diagnostic Check Tile (Green Tile) • Reflectance Black Glass • Sample Cup • Opaque Cover • Port Insert for Sample Cup • Port Plate - 31.8 mm (1.25 in) • Power Supply (90W) • Quick Start Guide • USB Flash Drive with Tomato License

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com